

Standard Serums

Adult Bovine Serum

Australian / New Zealand Origin

Bovogen Biologicals' Adult Bovine Serum is collected in a manner which minimises haemoglobin and endotoxin levels, and is processed in a validated sterile filtration facility.

Raw Material Origin

Bovine blood from adult bovine animals is collected from Australian and/or New Zealand registered export abattoirs under the supervision of Government Veterinary authorities.

In line with our quality systems, the origin of Bovogen ABS is traceable back to the abattoir, and is supported by the availability of full collection and manufacturing records.

Bovogen ABS origin is fully traceable, so we can guarantee our origin.

Appearance

Generally a clear, orange to brown solution.

Applications

Applications include:

- general additive in growth media for tissue and cell culture
- research or diagnostic assay reagent

Testing Information

GENERAL (Standard) VIRUSES / MYCOPLASMA (Standard)

(*) Bovine viruses only

Contact Bovogen for a Product Specification Sheet

Certificate of Suitability

Bovogen Biologicals Australian Origin / New Zealand ABS has a Certificate of Suitability (RO CEP 2014-142-Rev 00) from the European Directorate of Quality Medicines (EDQM).

Documentation

A Certificate of Analysis is supplied with all ABS orders. A Certificate of Origin is available in request.

Storage and Handling

Product is shipped and stored at -15 to -40°C.

Ordering Information

Contact Bovogen Biologicals regarding availability.